

Professional
Powder Equipment
Manufacturer



Powder
Equipment



Milling
Technology



Powder
Materials

TENCAN

Product Brochure



GLOVE BOX ACCESSORIES

Glove Box Dry Scroll Vacuum Pump

Oil-free dry scroll vacuum pump designed for glove box systems, providing clean and reliable vacuum performance without oil contamination for evacuation and controlled atmosphere applications.

<https://www.gloveboxs.com/products/glove-box-accessories/glove-box-dry-scroll-vacuum-pump.html>

Generated at 2026-08-17 04:04:51



Product Overview

Oil-free dry scroll vacuum pump designed for glove box systems, providing clean and reliable vacuum performance without oil contamination for evacuation and controlled atmosphere applications.

Technical Parameters

Parameter	Model 1
Parameter 1	
Parameter 2	



Product Introduction

Pumping speed □ 35 m³h⁻¹ / 21 ft³min⁻¹

Ultimate vacuum □ 1 x 10⁻² mbar / 8 x 10⁻³ Torr

Motor power □ 0.52kW

Inlet flange □ NW40

Electrical Requirements □ 115/230 V

Accessories & Customization

Accessories

Grinding jars, heating elements, sample holders, control modules and other matching accessories can be selected according to the product configuration.

Customization

For voltage, capacity, chamber size, process temperature or application requirements, please contact TENCAN for a suitable configuration.