

Professional
Powder Equipment
Manufacturer



Powder
Equipment



Milling
Technology



Powder
Materials

TENCAN

Product Brochure



GLOVE BOX ACCESSORIES

Glove Box Zirconia Oxygen Analyzer

Zirconia oxygen analyzer designed for continuous O₂ monitoring in glove box systems, providing stable and reliable measurement for controlled atmosphere laboratory and production applications.

<https://www.gloveboxs.com/products/glove-box-accessories/glove-box-zirconia-oxygen-analyzer.html>

Generated at 2026-08-17 04:05:44



Product Overview

Zirconia oxygen analyzer designed for continuous O₂ monitoring in glove box systems, providing stable and reliable measurement for controlled atmosphere laboratory and production applications.

Technical Parameters

Parameter	Model 1
Parameter 1	
Parameter 2	



Product Introduction

Range: 0-1000 ppm

Power Supply: DC 24V

Signal: 4-20 mA

Accuracy: Within 100 ppm, ± 1 ppm or 2% of the displayed value

Display: Via PLC or independent display

Accessories & Customization

Accessories

Grinding jars, heating elements, sample holders, control modules and other matching accessories can be selected according to the product configuration.

Customization

For voltage, capacity, chamber size, process temperature or application requirements, please contact TENCAN for a suitable configuration.