

Professional
Powder Equipment
Manufacturer



Powder
Equipment



Milling
Technology



Powder
Materials

TENCAN

Product Brochure



GLOVE BOX ACCESSORIES

Gas and Liquid Feedthrough

Gas and liquid feedthrough designed for connecting equipment inside a glove box or vacuum enclosure to external gas and liquid lines. Features an 8 mm stainless steel tube and KF flange sealing for reliable vacuum service.

<https://www.gloveboxs.com/products/glove-box-accessories/gas-liquid-feedthrough.html>

Generated at 2026-08-17 04:07:55



Product Overview

Gas and liquid feedthrough designed for connecting equipment inside a glove box or vacuum enclosure to external gas and liquid lines. Features an 8 mm stainless steel tube and KF flange sealing for reliable vacuum service.

Technical Parameters

Parameter	Model 1
Parameter 1	
Parameter 2	



Product Introduction

Medium: Gas, Liquid

Pipeline: $\phi 8$ Stainless Steel Tubing

Seal: KF Flange Seal

Leak Rate: 10-Pa·L/s

Application: Mainly used for connecting equipment inside the glove box vacuum chamber to external equipment.

Accessories & Customization

Accessories

Grinding jars, heating elements, sample holders, control modules and other matching accessories can be selected according to the product configuration.

Customization

For voltage, capacity, chamber size, process temperature or application requirements, please contact TENCAN for a suitable configuration.