

Professional
Powder Equipment
Manufacturer



Powder
Equipment



Milling
Technology



Powder
Materials

TENCAN

Product Brochure



GLOVE BOX ACCESSORIES

Glove Box Butyl Rubber Gloves

Butyl rubber gloves designed for glove box systems, providing excellent airtightness, chemical resistance and reliable protection for controlled atmosphere applications.

<https://www.gloveboxs.com/products/glove-box-accessories/glove-box-butyl-rubber-gloves.html>

Generated at 2026-08-17 04:10:35



Product Overview

Butyl rubber gloves designed for glove box systems, providing excellent airtightness, chemical resistance and reliable protection for controlled atmosphere applications.

Technical Parameters

Parameter	Model 1
Parameter 1	
Parameter 2	



Product Introduction

Caliber: 8 inches/220mm

Thickness: 0.4mm/0.8mm

Length: 800mm

Material: Butyl

Accessories & Customization

Accessories

Grinding jars, heating elements, sample holders, control modules and other matching accessories can be selected according to the product configuration.

Customization

For voltage, capacity, chamber size, process temperature or application requirements, please contact TENCAN for a suitable configuration.